

Using IDEXX VetLab with Cornerstone Software


This guide explains the recommended way to create an IDEXX VetLab* in-clinic laboratory request and view the test results.

IMPORTANT: This document applies to Cornerstone* Software 8.4 NEXT or later. Some features may have changed from previous versions. Refer to additional online resources and Cornerstone help [?](#) for more information.

Create a laboratory request using the Patient Visit List or an invoice

Note: The invoice item must be linked to the profile, and the **Use special actions** check box must be selected in the INCLINIC laboratory preferences. For details, see the [IDEXX VetLab with Cornerstone Software Setup Guide](#).

1. On the **Patient Clipboard*** , in the **Patient List** area, right-click the patient name, and then select either **Patient Visit List (PVL)** or **Invoice**.
2. Enter the invoice item ID for the test, and then press **TAB** to open a new laboratory request.



Item ID	Item Description	Profile/test	Price	Center
1	7200 CATALYST CHEM 17 CLIP	Chem17 CLIP	\$79.33	H
2				
3				

3. Enter or modify the staff information as needed.
4. Enter values in the **Drawn by** box, if required, and click **OK**.
Note: If you have multiple IDEXX VetLab Stations connected to the Cornerstone software, select the station you want to use.
5. If more tests are needed, click **OK**; at the Print Box, click **Close** and then add the additional invoice items following **steps 2-4**.
6. After all tests are entered, in the last Print Box, update the number of label copies as needed, and then click **Print**.

Note: The steps above may vary, depending on your laboratory default settings.

Running tests on the IDEXX VetLab Station


For detailed information about using the IDEXX VetLab Station, see the [IDEXX VetLab Station Quick Reference Guide](#).

1. Tap the desired patient from the Pending or **Census** list on the Home screen and tap **Run**.
2. Confirm/specify the patient/practice information, select the desired instrument icon(s), and tap **Run**.

View results in the patient record using VetConnect* PLUS

If VetConnect* PLUS has been activated, the VetConnect PLUS tab will be available whenever you view laboratory results.

1. On the toolbar, click **Patient Clipboard** . VetConnect* PLUS results can also be accessed from the Diagnostics Tab of the Daily Planner.
2. Select the client and patient record.
3. Select the **VetConnect PLUS** tab. For more information on using VetConnect PLUS with Cornerstone, press **F1**.



Note: An internet connection is required to view VetConnect PLUS results. VetConnect PLUS display and functionality will vary on versions earlier than 9.3.

Additional guide

[Laboratory daily checklist for Cornerstone software](#)

Find this guide

Look under “Resources by Topic” on the [Cornerstone Help](#) page.