SNAP® Feline proBNP

In vitro test for the measurement of circulating NTproBNP in feline serum and EDTA plasma.

Precautions and Warnings

- All wastes should be properly decontaminated prior to disposal.
- Do not mix components from kits with different serial numbers.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- WARNING:** Conjugate – H316/P332+P313/EUH208. Causes mild skin irritation. If skin irritation occurs: Get medical advice/attention. Contains Kathon. May produce an allergic reaction.

Storage

- Store at 2–8°C until expiration date.

Kit Components

Item	Reagents	Quantity
1	1 bottle of anti-chicken:HRPO/anti-NTproBNP: HRPO conjugate (Green)	3,0 mL
2	SNAP Device	5 or 10
<i>Each SNAP device contains 0.4 mL of wash solution and 0.6 mL of substrate solution.</i>		
Other Components: transfer pipettes, sample tubes and reagent rack		

Sample Information

- Samples must be at room temperature (18–25°C) before beginning the test procedure.
- Serum or EDTA plasma, either fresh or stored at 2–8°C for up to one week, can be used.
- For longer storage, serum and EDTA plasma can be frozen (-20°C or colder) and then centrifuged before use.
- Hemolyzed, icteric or lipemic samples will not affect test results.

Test Procedure

- Allow all components to equilibrate at room temperature (18–25°C) for at least 30 minutes. **Do not heat.**

- Hold the required transfer pipette vertically and dispense **3 drops of sample** into a new sample tube.
- Holding the bottle vertically, add **5 drops of green conjugate** to the sample tube.
- Cap the sample tube and mix thoroughly by inverting it 3–5 times.
- Place the device on a horizontal surface. Add the entire contents of the sample tube to the sample well, being careful not to splash the contents outside of the sample well.


The sample will flow across the result window, reaching the activation circle in 30–60 seconds. Some sample may remain in the sample well.

- When color **FIRST** appears in the activation circle, push the activator firmly until it is flush with the device body.

Note: Some samples may not flow to the activation circle within 60 seconds. In this case, press the activator after the sample has flowed across the result window.

- Read test result **10 minutes** from the time of activation.

Interpreting the Test Results

Abnormal Level

If the color intensity of the sample spot is equal to or darker than the color intensity of the reference spot, NTproBNP levels are abnormal.

Normal Level

If the color intensity of the sample spot is lighter than the color intensity of the reference spot, NTproBNP levels are normal.

Note: It is possible for no color to develop on the sample spot.

Invalid Results

- Background**—If the sample is allowed to flow past the activation circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No Color Development**—If reference spot does not develop color, repeat the test.

Percent agreement to Cardiopet

402 samples were tested on both SNAP Feline proBNP and Cardiopet proBNP tests. Visual results from SNAP Feline proBNP demonstrated 99.5% (95% CI: 97.1 – 99.9) agreement with Cardiopet proBNP if the sample had a concentration of NT-proBNP < 100 pmol/L (n=187) and 95.0% (95% CI: 90.0 – 97.5) agreement with Cardiopet proBNP if the sample concentration of NT-proBNP was ≥ 270 pmol/L (n=140).


IDEXX SNAPshot Dx* Analyzer

Test data can also be read using the SNAPshot Dx* Analyzer. A complete description of how to enter patient data and read test results using the SNAPshot Dx Analyzer can be found in the SNAPshot Dx Analyzer user guide.

Symbol Descriptions

	Use by date		Authorized Representative in the European Community
	Batch Code (Lot)		Consult instructions for use
	Serial number		In vitro diagnostic
	Temperature limitation		Date of manufacture
	Manufacturer		Catalog Number

IDEXX Customer Support
USA/Canada 1 800 248 2483 • Europe 00800 1234 3399 • Australia 1300 44 33 99
idexx.com

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* SNAP is a trademark or registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.
© 2015 IDEXX Laboratories, Inc. All rights reserved. 06-0001179-01
© 2015 IDEXX Laboratories, Inc. Alle rechten voorbehouden.
IDEXX Europe B.V.