


SNAP® T₄ Testing Guide

Canine



SNAP® T₄ Results

| | | |
|----------------|---------------|----------------|
| Low | <0.8 µg/dL | (<10 nmol/L) |
| Borderline Low | 0.8–1.5 µg/dL | (10–20 nmol/L) |
| Normal | 1.6–5.0 µg/dL | (21–64 nmol/L) |
| High | >5.0 µg/dL | (>64 nmol/L) |
| Therapeutic | 3.0–6.0 µg/dL | (39–77 nmol/L) |

Canine Screening

- Dogs with subnormal T₄ results may be hypothyroid or euthyroid sick.
- Dogs with T₄ results in the borderline low range may be hypothyroid.
- In dogs with consistent clinical signs, consider a free T₄ and endogenous TSH to aid in confirming hypothyroidism.


Hypothyroidism Therapeutic Monitoring

For dogs on thyroid supplement, acceptable 4–6-hour post-pill total T₄ values will generally fall in the upper end of the normal range or slightly above.

Note: 1 µg/dL is equal to 12.87 nmol/L. A result that falls within the borderline range of the assay should be considered ambiguous.

SNAP[®] T₄ Testing Guide

Feline



SNAP[®] T₄ Results

| | | |
|-----------------|---------------|----------------|
| Low | <1.0 µg/dL | (<13 nmol/L) |
| Normal | 1.0–5.0 µg/dL | (13–64 nmol/L) |
| Borderline High | 2.5–5.0 µg/dL | (32–64 nmol/L) |
| High | >5.0 µg/dL | (>64 nmol/L) |

Feline Screening

Older cats with consistent clinical signs and T₄ values in the borderline high range may have early hyperthyroidism or a concurrent nonthyroidal illness. In these cases, consider a free T₄, a T₃ suppression test or radionuclide thyroid imaging to aid in confirming the diagnosis.


Hypothyroidism Therapeutic Monitoring

Following treatment with methimazole, T₄ values will generally fall within the lower end of the normal range.

Note: 1 µg/dL is equal to 12.87 nmol/L. A result that falls within the borderline range of the assay should be considered ambiguous.

SNAP[®] T₄ Testing Guide

Equine



SNAP[®] T₄ Results

| | | |
|--------|---------------|----------------|
| Low | <0.9 µg/dL | (<12 nmol/L) |
| Normal | 0.9–2.8 µg/dL | (12–36 nmol/L) |
| High | >2.8 µg/dL | (>36 nmol/L) |

Note: 1 µg/dL is equal to 12.87 nmol/L.