
# Enterolert\*-DW

06-18085-07

**IDEXX**

IDX 33/03 - 10/13  
WATER ANALYSIS METHODS  
<http://nf-validation.afnor.org>

The method Enterolert\*-DW/Quanti-Tray\* for water analysis is granted NF Validation by AFNOR Certification as an alternative method to the standard ISO 7899-2 for enumeration of intestinal enterococci in drinking water, under the Certificate number: IDX 33/03 – 10/13. For more information about end of validity, please refer to the certificate NF Validation available on website mentioned above.

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IDEXX Water Quality Control Laboratory is accredited to ISO/IEC 17025:2017

**IDEXX**

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA  
[idexx.com/water](http://idexx.com/water)

**For Technical Support, please call:**

North/South America: 1 207 556 4496/1 800 321 0207  
Europe: 00800 4339 9111  
UK: +44 (0) 1638 676800  
China: +86 21 61279528  
Japan: 03 5301 6800  
Australia: 1300 443 399

# Enterolert\*-DW Test Kit

## Introduction

Enterolert\*-DW is designed to detect *Enterococcus* species in drinking water samples. It is based on IDEXX's proprietary Defined Substrate Technology\* (DST\*). Enterolert-DW, coupled with the IDEXX Quanti-Tray\* System, provides quantitative confirmed results in 24 hours. Enterolert-DW utilizes Ortho-Nitrophenyl- $\beta$ -D-glucoside as a nutrient indicator and incorporates a specifically designed blue background colour in its formulation. When the substrate is metabolized by enterococci, the sample turns from blue to green to indicate a positive detection. Any change from the original colour to green is considered a positive result. No ultraviolet light source is required. Enterolert-DW detects enterococci in drinking water samples in 24 hours.

## Storage

Store at 2–25°C away from direct light and humidity.

## Quanti-Tray\* Enumeration Procedure

1. Add contents of one pack to a 100 mL water sample in a sterile vessel.
2. Cap vessel and shake until dissolved.
3. Pour sample/reagent mixture into a Quanti-Tray\* or Quanti-Tray\*/2000 and seal with an IDEXX Quanti-Tray\* Sealer.
4. Place the sealed tray in a 41 ± 0.5°C incubator for 24 hours.
5. Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

## Result Interpretation

Appearance	Result
Blue	Negative for enterococci
Green	Positive for enterococci

- Enterolert-DW results are definitive at 24–28 hours. In addition, positive wells for enterococci observed before 24 hours and negative wells observed after 28 hours are also valid.

## Procedural Notes

- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, a negative control can be used when interpreting results.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- Enterolert-DW can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*<sup>1</sup> MPN tables should be used to find Most Probable Numbers.
- Enterolert-DW is a primary drinking water test. Enterolert-DW performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- Aseptic technique should always be followed when using Enterolert-DW. Dispose of in accordance with Good Laboratory Practice.

## Quality Control Procedures

1. One of the following quality control procedures is recommended for each lot of Enterolert-DW:
  - A. IDEXX-QC Enterococci DW<sup>2</sup>: *Enterococcus faecalis*, *Serratia marcescens*.
  - B.
    - i. For each of the American Type Culture Collection (ATCC)<sup>3</sup> bacterial strains (*Enterococcus faecium* 35667 and *Serratia marcescens* 43862), streak the culture onto labeled TSA or Blood Agar plates and incubate at 35 ± 2°C for 18–24 hours.
    - ii. For positive controls, prepare proper dilutions to reach a concentration of approximately 60–80 CFU/100 mL.
    - iii. For non-target strains, prepare proper dilutions to a suggested inoculation up to 10<sup>6</sup> CFU/100 mL.
2. Results should match the Result Interpretation table above.

**NOTE:** IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at [idexx.com/water](http://idexx.com/water).

1. Eaton AD, Clesceri LS, Rice EW, Greenberg AE, Franson MAH, eds. *Standard Methods for the Examination of Water and Wastewater*. 21st ed. Washington, DC: American Public Health Association; 2005.

2. IDEXX-QC Enterococci-DW—IDEXX Catalog #UN3373-WQC-ENT DW

3. American Type Culture Collection 1-800-638-6597 [atcc.org](http://atcc.org)

\*Enterolert, Defined Substrate Technology, DST and Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Patent information: [idexx.com/patents](http://idexx.com/patents).

© 2019 IDEXX Laboratories, Inc. All rights reserved.

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© 2019 IDEXX Laboratories, Inc. Wszelkie prawa zastrzeżone.