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Références

- Keene BW, Atkins CE, Bonagura JD, *et al.* ACVIM consensus guidelines for the diagnosis and treatment of myxomatous mitral valve disease in dogs. *J Vet Intern Med.* 2019;33(3):1127-1140. doi:10.1111/jvim.15488
- Luis Fuentes V, Abbott J, Chetboul V, *et al.* ACVIM consensus statement guidelines for the classification, diagnosis, and management of cardiomyopathies in cats. *J Vet Intern Med.* 2020;34(3):1062-1077. doi:10.1111/jvim.15745
- Bode L, Lau S. Interpreting cardiac biomarkers in small animal practice. *In Pract.* 2023;45(7):372-382. doi:10.1002/inpr.349
- Boswood A. Biomarkers in cardiovascular disease: beyond natriuretic peptides. *J Vet Cardiol.* 2009;11 Suppl 1:S23-S32. doi:10.1016/j.jvc.2009.01.003
- Cahill RJ, Pigeon K, Strong-Townsend MI, Drexel JP, Clark GH, Buch JS. Analytical validation of a second-generation immunoassay for the quantification of N-terminal pro-B-type natriuretic peptide in canine blood. *J Vet Diagn Invest.* 2015;27(1):61-67. doi:10.1177/1040638714562826
- Mainville CA, Clark GH, Esty KJ, *et al.* Analytical validation of an immunoassay for the quantification of N-terminal pro-B-type natriuretic peptide in feline blood. *J Vet Diagn Invest.* 2015;27(4):414-421. doi:10.1177/1040638715588330
- Machen MC, Oyama MA, Gordon SG, *et al.* Multi-centered investigation of a point-of-care NT-proBNP ELISA assay to detect moderate to severe occult (pre-clinical) feline heart disease in cats referred for cardiac evaluation. *J Vet Cardiol.* 2014;16(4):245-255. doi:10.1016/j.jvc.2014.09.002
- Hezzell MJ, Rush JE, Humm K, *et al.* Accuracy of history, physical examination, cardiac biomarkers, and biochemical variables in identifying dogs with stage B2 degenerative mitral valve disease. *J Vet Intern Med.* 2016;30(2):536-542. doi:10.1111/jvim.13831
- Lu TL, Côté E, Kuo YW, Wu HH, Wang WY, Hung YW. Point-of-care N-terminal pro B-type natriuretic peptide assay to screen apparently healthy cats for cardiac disease in general practice. *J Vet Intern Med.* 2021;35(4):1663-1672. doi:10.1111/jvim.16156
- Carter R, Partington C, Novo Matos J. Accuracy of combined signalment, physical examination and cardiac biomarkers in screening cats with asymptomatic murmurs for cardiac disease. *J Vet Intern Med.* 2026;40(1):aala037. doi:10.1093/jvims/aala037
- Wilshaw J, Rosenthal SL, Wess G, *et al.* Accuracy of history, physical examination, cardiac biomarkers, and biochemical variables in identifying dogs with stage B2 degenerative mitral valve disease. *J Vet Intern Med.* 2021;35(2):755-770. doi:10.1111/jvim.16083
- Wesselowski S, Gordon SG, Fries R, *et al.* Use of physical examination, electrocardiography, radiography, and biomarkers to predict echocardiographic stage B2 myxomatous mitral valve disease in preclinical Cavalier King Charles Spaniels. *J Vet Cardiol.* 2023;50:1-16. doi:10.1016/j.jvc.2023.10.001
- Fox PR, Oyama MA, Hezzell MJ, *et al.* Relationship of plasma N-terminal pro-brain natriuretic peptide concentrations to heart failure classification and cause of respiratory distress in dogs using a 2nd generation ELISA assay. *J Vet Intern Med.* 2015;29(1):171-179. doi:10.1111/jvim.12472
- Oyama MA, Rush JE, Rozanski EA, *et al.* Assessment of serum N-terminal pro-B-type natriuretic peptide concentration for differentiation of congestive heart failure from primary respiratory tract disease as the cause of respiratory signs in dogs. *JAVMA.* 2009;235(11):1319-1325. doi:10.2460/javma.235.11.1319
- Fox PR, Oyama MA, Reynolds C, *et al.* Utility of plasma N-terminal pro-brain natriuretic peptide (NT-proBNP) to distinguish between congestive heart failure and non-cardiac causes of acute dyspnea in cats. *J Vet Cardiol.* 2009;11 Suppl 1:S51-S61. doi:10.1016/j.jvc.2008.12.001
- Humm K, Hezzell M, Sargent J, Connolly DJ, Boswood A. Differentiating between feline pleural effusions of cardiac and non-cardiac origin using pleural fluid NT-proBNP concentrations. *J Small Anim Pract.* 2013;54(12):656-661. doi:10.1111/jsap.12152
- Singletary GE, Rush JE, Fox PR, Stepien RL, Oyama MA. Effect of NT-pro-BNP assay on accuracy and confidence of general practitioners in diagnosing heart failure or respiratory disease in cats with respiratory signs. *J Vet Intern Med.* 2012;26(3):542-546. doi:10.1111/j.1939-1676.2012.00916.x
- Dukes-McEwan J, Garven KE, Lopez Alvarez J, Oliveira P, Molskula PF, Willis R. Usefulness of cardiac biomarker screening to detect dilated cardiomyopathy in Doberman Pinschers. *J Small Anim Pract.* 2022;63(4):275-285. doi:10.1111/jsap.13455
- Morey AGN, Lamb KE, Karnia JJ, *et al.* N-terminal brain natriuretic peptide, cardiac troponin-I, and point-of-care ultrasound in dogs with cardiac and noncardiac causes of nonhemorrhagic ascites. *J Vet Intern Med.* 2023;37(3):900-909. doi:10.1111/jvim.16702
- Borgarelli M, Ferasin L, Lamb K, *et al.* The predictive value of clinical, radiographic, echocardiographic variables and cardiac biomarkers for assessing risk of the onset of heart failure or cardiac death in dogs with preclinical myxomatous mitral valve disease enrolled in the DELAY study. *J Vet Cardiol.* 2021;36:77-88. doi:10.1016/j.jvc.2021.04.009
- Kaplan JL, Visser LC, Gunther-Harrington CT, *et al.* Effect of standard-dose and high-dose pimobendan on select indices of renal and cardiac function in dogs with American College of Veterinary Internal Medicine stage B2 myxomatous mitral valve disease. *J Vet Intern Med.* 2022;36(6):1892-1899. doi:10.1111/jvim.16537
- Hezzell MJ, Block CL, Laughlin DS, Oyama MA. Effect of prespecified therapy escalation on plasma NT-proBNP concentrations in dogs with stable congestive heart failure due to myxomatous mitral valve disease. *J Vet Intern Med.* 2018;32(5):1509-1516. doi:10.1111/jvim.15228
- Wess G, Daisenberger P, Mahling M, Hirschberger J, Hartmann K. Utility of measuring plasma N-terminal pro-brain natriuretic peptide in detecting hypertrophic cardiomyopathy and differentiating grades of severity in cats. *Vet Clin Pathol.* 2011;40(2):237-244. doi:10.1111/j.1939-165X.2011.00305.x
- Pierce KV, Rush JE, Freeman LM, Cunningham SM, Yang VK. Association between survival time and changes in NT-proBNP in cats treated for congestive heart failure. *J Vet Intern Med.* 2017;31(3):678-684. doi:10.1111/jvim.14690