SNAP FeLV Antigen Test

Test accuracy

Sensitivity and specificity of the SNAP FeLV Antigen Test

Analyte	Reference method	SNAP test result		Total	Sensitivity (95% CL)
		+	–		Specificity (95% CL)
FeLV[1]	+	158	0	407	100% (97.7%–100%)
	–	2	247		99.2% (97.1-99.9%)

Source: Singh S, Davenport KA, Schooley E, Ruggiero A, Nassar S, Buch J, Chandrashekar R. Diagnostic Accuracy of a Point-of-Care Immunoassay for Feline Immunodeficiency Virus Antibodies, Feline Leukemia Virus Antigen, and Dirofilaria immitis Antigen. Viruses. 2023 Oct 19;15(10):2117. doi:10.3390/v15102117. PMID:37896893; PMCID: PMC106112086.

Reference method

1. PetChek* Feline Leukemia Virus Antigen Test

For interpretation and diagnostic algorithms, visit **idexx.com/FeLVAntigenTestResources**.

© 2026 IDEXX Laboratories, Inc. All rights reserved. • 09-2693053-00 • The IDEXX Privacy Policy is available at idexx.com.

All trademarks are owned by IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.